Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

MK-571 sodium salt

Product overview

Name MK-571 sodium salt

 Cat No
 HB3976

 Purity
 >97%

Description Selective leukotriene LTD4 receptor antagonist. Potent MRP1 inhibitor.

Biological Data

Biological description Selective leukotriene LTD4 receptor antagonist. Potent MRP1 inhibitor. Antinociceptive.

Solubility & Handling

Solubility overview Soluble in water (25mM) and in EtOH (10mM)

Storage instructions +4°C

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

Chemical Data

Molecular Weight 537

Chemical structure

InChiKey XNAYQOBPAXEYLI-AAGWESIMSA-M

Appearance Light yellow solid