Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



# **DATASHEET**

S14161

#### **Product overview**

 Name
 \$14161

 Cat No
 HB4068

 Purity
 >99%

**Description** pan Class I PI3K isoforms inhibitor

### **Biological Data**

Biological description pan Class I PI3K isoforms inhibitor. Potent antileukemia and antimyeloma compound. Apoptosis

inducer. Delays tumor growth. Cyclin D2 promoter transactivation inhibitor.

## **Solubility & Handling**

Storage instructions +4°C

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

#### **Chemical Data**

Molecular Weight Chemical structure 315.3 CH<sub>3</sub> NO<sub>2</sub>

**CAS Number** 883046-50-2

InChiKey FZXDHIUPNWOMT-UHFFFAOYSA-N

Appearance Yellow solid