

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

RIPA buffer (10x)

Product overview

Name	RIPA buffer (10x)
Cat No	HB6569
Biological description	10X RIPA lysis buffer for cell lysis and protein solubilization for both adherent and suspension cultured mammalian cells. Effectively extracts cytoplasmic, nuclear and membrane proteins and is suitable for applications including SDS-PAGE, Western blot, immunoprecipitation, ELISA and BCA assays. Preparation: Dilute to 1X in deionized water
Biological action Description	Buffer 10X RIPA lysis buffer

Biological Data

Application notes	<ul style="list-style-type: none">• Please note this is a 10x buffer. To dilute to 1x add 100ml of RIPA buffer (10x) to 900ml of dH₂O• Please see our guide on how to prepare samples for Western blots for using this RIPA buffer to prepare protein samples for use in SDS-PAGE.
--------------------------	--

Solubility & Handling

Storage instructions Important	Room temperature This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.
---------------------------------------	---
