Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



# **DATASHEET**

TBS-T with Tween 20 (20x)

#### **Product overview**

Name TBS-T with Tween 20 (20x)

Cat No HB6971

**Biological description** 20x Tris Buffered Saline wash buffer for use during immunolabeling techniques such as IHC. When

diluted to 1x contains 130mM NaCl, 20mM Tris, 0.1% Tween 20, 0.002% ProClin-300, pH7.4.

Biological action Buffer

**Description** 20x Tris Buffered Saline wash buffer to remove non-specific binding and for use during

immunolabeling techniques such as IHC.

## **Biological Data**

Application notes This is a 20x buffer: to use, dilute in dH<sub>2</sub>O (50ml of 20x buffer + 950ml dH<sub>2</sub>O)

### **Solubility & Handling**

Storage instructions Room temperature

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

#### **Chemical Data**

Appearance Colorless solution