

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

TBS-T with Tween 20 (20x)

Product overview

Name	TBS-T with Tween 20 (20x)
Cat No	HB6971
Biological description	20x Tris Buffered Saline wash buffer for use during immunolabeling techniques such as IHC. When diluted to 1x contains 130mM NaCl, 20mM Tris, 0.1% Tween 20, 0.002% ProClin-300, pH7.4.
Biological action	Buffer
Description	20x Tris Buffered Saline wash buffer to remove non-specific binding and for use during immunolabeling techniques such as IHC.

Biological Data

Application notes	This is a 20x buffer: to use, dilute in dH ₂ O (50ml of 20x buffer + 950ml dH ₂ O)
--------------------------	--

Solubility & Handling

Storage instructions	Room temperature
Important	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

Chemical Data

Appearance	Colorless solution
-------------------	--------------------
