

Hello Bio, Inc.  
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500  
F. 609-228-4994

customercare-usa@hellobio.com



---

## DATASHEET

### TBS-T with Tween 20 (20x)

---

#### Product overview

<b>Name</b>	TBS-T with Tween 20 (20x)
<b>Cat No</b>	HB6971
<b>Biological description</b>	20x Tris Buffered Saline wash buffer for use during immunolabeling techniques such as IHC. When diluted to 1x contains 130mM NaCl, 20mM Tris, 0.1% Tween 20, 0.002% ProClin-300, pH7.4.
<b>Biological action</b>	Buffer
<b>Description</b>	20x Tris Buffered Saline wash buffer to remove non-specific binding and for use during immunolabeling techniques such as IHC.

---

#### Biological Data

<b>Application notes</b>	This is a 20x buffer: to use, dilute in dH <sub>2</sub> O (50ml of 20x buffer + 950ml dH <sub>2</sub> O)
--------------------------	--

---

#### Solubility & Handling

<b>Storage instructions</b>	Room temperature
<b>Important</b>	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

---

#### Chemical Data

<b>Appearance</b>	Colorless solution
-------------------	--------------------

---