

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

Product overview

| | |
|-------------------------------|--|
| Name | Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0) |
| Cat No | HB7943 |
| Biological description | Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue sections. Optimal for use with primary antibodies that require citrate buffer (pH 6.0) pretreatment. This product contains detergent for emulsification of the paraffin. 1X Dilution: The 100X stock solution should be diluted 100-fold with distilled water before use. |
| Biological action | Buffer |
| Applications | IHC, IHC-P |
| Description | Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue sections |

Solubility & Handling

| | |
|-----------------------------|---|
| Storage instructions | Room temperature |
| Important | This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use. |
