Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

Product overview

Name Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

Cat No HB7943

Biological description Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue

sections. Optimal for use with primary antibodies that require citrate buffer (pH 6.0) pretreatment. This product contains detergent for emulsification of the paraffin. 1X Dilution: The 100X stock solution

should be diluted 100-fold with distilled water before use.

Biological action ApplicationsBuffer
IHC, IHC-P

Description Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue

sections

Solubility & Handling

Storage instructions Room temperature

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.