

Hello Bio, Inc.  
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500  
F. 609-228-4994

customercare-usa@hellobio.com



---

## DATASHEET

Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)

---

### Product overview

|                               |  |
|-------------------------------|--|
| <b>Name</b>                   | Antigen Retrieval Buffer (100X Tris-EDTA Buffer, pH 9.0)   |
| <b>Cat No</b>                 | HB7943   |
| <b>Biological description</b> | Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue sections. Optimal for use with primary antibodies that require citrate buffer (pH 6.0) pretreatment. This product contains detergent for emulsification of the paraffin. 1X Dilution: The 100X stock solution should be diluted 100-fold with distilled water before use. |
| <b>Biological action</b>      | Buffer   |
| <b>Applications</b>           | IHC, IHC-P   |
| <b>Description</b>            | Enables one-step rehydration and antigen retrieval in formalin-fixed, paraffin-embedded tissue sections  |

---

### Solubility & Handling

|                             |   |
|-----------------------------|---|
| <b>Storage instructions</b> | Room temperature  |
| <b>Important</b>            | This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use. |

---