Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

PBS buffer with Tween 20 (20x)

Product overview

Name PBS buffer with Tween 20 (20x)

Cat No HB8088

Biological description 20x wash buffer for use during immunolabeling techniques such as IHC. Tween20 promotes more

effective washing to decrease non-specific background staining. When diluted to 1x contains 0.1%

Tween-20 in PBS, pH7.4.

Species of origin Synthetic Biological action Buffer

Description 20x wash buffer for use during immunolabeling techniques such as IHC.

Solubility & Handling

Storage instructions Room temperature

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.