

Hello Bio, Inc.  
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500  
F. 609-228-4994

customercare-usa@hellobio.com



---

## DATASHEET

### PBS buffer with Tween 20 (20x)

---

#### Product overview

<b>Name</b>	PBS buffer with Tween 20 (20x)
<b>Cat No</b>	HB8088
<b>Biological description</b>	20x wash buffer for use during immunolabeling techniques such as IHC. Tween20 promotes more effective washing to decrease non-specific background staining. When diluted to 1x contains 0.1% Tween-20 in PBS, pH7.4.
<b>Species of origin</b>	Synthetic
<b>Biological action</b>	Buffer
<b>Description</b>	20x wash buffer for use during immunolabeling techniques such as IHC.

---

#### Solubility & Handling

<b>Storage instructions</b>	Room temperature
<b>Important</b>	This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.

---