

Hello Bio, Inc.
304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500
F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

PBS buffer with Tween 20 (20x)

Product overview

| | |
|-------------------------------|---|
| Name | PBS buffer with Tween 20 (20x) |
| Cat No | HB8088 |
| Biological description | 20x wash buffer for use during immunolabeling techniques such as IHC. Tween20 promotes more effective washing to decrease non-specific background staining. |
| Biological action | Buffer |
| Description | 20x wash buffer for use during immunolabeling techniques such as IHC. |

Solubility & Handling

| | |
|-----------------------------|---|
| Storage instructions | Room temperature |
| Important | This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use. |
