Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

Antibody Diluent

Product overview

Name Antibody Diluent

Cat No HB8223

Biological description Diluent used to dilute primary antibodies to suitable working concentrations (e.g. for IHC) or for

storage. Reduces background and non-specific binding of the primary antibody.

Biological action Buffer

Applications ELISA, ICC, IHC, WB, IHC(IF), Histoblot

Description Diluent used to dilute primary antibodies or for storage

Solubility & Handling

Storage instructions Room temperature

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.

Chemical Data

Formulation Liquid