Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

Citrate Buffer (10x) pH 6.0

Product overview

Name Citrate Buffer (10x) pH 6.0

Cat No HB8687

Biological descriptionCitrate buffer (pH6) for antigen retrieval. Used for antigen retrieval prior to IHC (immunohistochemistry)

on formalin-fixed, paraffin-embedded tissue sections mounted on glass slides to increase staining

intensity with many primary antibodies.

Description Citrate buffer (pH6) for target retrieval

Solubility & Handling

Storage instructions Room temperature

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.