Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

RIPA buffer

Product overview

Name RIPA buffer Cat No HB9619

Biological description Ready to use cell and tissue lysis buffer suitable for lysing cells for radioIP. Also results in low

background in immunoprecipitation/pull down assays

Biological action Buffer

Description Ready to use cell and tissue lysis buffer

Biological Data

Application notes Please see our guide on how to prepare samples for Western blots for using this RIPA buffer to

prepare protein samples for use in SDS-PAGE.

Solubility & Handling

Storage instructions Room temperature

Storage of solutions Prepare and use solutions on the same day if possible. Store solutions at -20 °C for up to one month if

storage is required. Equilibrate to RT and ensure the solution is precipitate free before use.

Shipping Conditions Stable for ambient temperature shipping. Follow storage instructions on receipt.

Important This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not

for human or veterinary use.