Hello Bio, Inc. 304 Wall St., Princeton, NJ 08540 USA

T. 609-683-7500 F. 609-228-4994

customercare-usa@hellobio.com



DATASHEET

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ValidAb™

Product overview

Name Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ValidAb™

 Cat No
 HB9524

 Host
 Goat

 Clonality
 Polyclonal


 Target
 Rabbit IgG H&L

Conjugate Horseradish peroxidase (HRP)

Description Goat Anti-Rabbit IgG H&L Horseradish peroxidase (HRP) secondary antibody. Part of the ValidAb

TM range of highly validated, data-rich antibodies.

Validation data



Product information

Immunogen Purified rabbit IgG

Isotype

Purification notes Immunogen affinity chromatography. Pre-adsorbed with mouse and human serum proteins

Concentration 0.35mg/ml

Formulation 50% Glycerol, 50% PBS pH 7.4 with 0.05% ProClin-300

Tested applications

Applications ELISA, IHC, WB

Western blot optimal Dependent on assay. A 1:20,000 dilution is a good starting point. concentration

Open data link Please follow this link to OSF

Storage & Handling

Storage instructions -20 °C Shipping Conditions On ice

References

Coordinated interplay between palmitoylation, phosphorylation and SUMOylation regulates kainate receptor surface expression.

Yucel BP et al (2023) Frontiers in molecular neuroscience 16

PubMedID 37868810

The induction of p53 correlates with defects in the production, but not the levels, of the small ribosomal subunit and stalled large ribosomal subunit biogenesis.

Eastham MJ et al (2023) Nucleic acids research 51 **PubMedID** 37526268

Western blot: technique, theory, and trouble shooting.

Mahmood T et al (2012) North American journal of medical sciences 4

PubMedID 23050259