

Certificate of Analysis



Product Identification

Catalog number	HB6565
Compound name	Recombinant mouse MANF protein
Description	Produced in E. Coli as a single, non-glycosylated, polypeptide chain of 158 Aas. Molecular weight of 18.2 kDa
Batch number	E2213-1-1
Batch molecular weight	18.2 kDa

Method	QC requirement	QC Result
HPLC	Reverse Phase HPLC shows >98% purity	Meets specification: 98.0%
Physical appearance	White lyophilised powder solid	Meets specification: White solid
Solubility	Soluble in aqueous buffer	Meets specification

Produced by

Signature

Passed by

Signature

Date

22/06/2023

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.