

Certificate of Analysis



Product Identification

Catalog number HB8373
Compound name Recombinant human NT-3 protein
Batch number E1823-1-1
Batch molecular weight 27.2 kDa

Method	QC requirement	QC Result
HPLC Peptide content	Reverse Phase HPLC shows >95% purity Peptide content between 60 % and 100 %	Meets specification: 98.0% Meets specification: 100%
Physical appearance Solubility	White to off-white solid Soluble in water (0.1 mg/ml)	Meets specification: White solid Meets specification

Produced by Richard Patterson

Signature

Passed by Steve Roome

Signature

Date

09-12-21

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.