

Certificate of Analysis



Product Identification

Catalog number	HB9314
Compound name	Recombinant human Persephin / PSPN protein
Description	PSPN human recombinant produced in E. Coli is a disulphide-linked homodimer containing 2x96 AA and a molecular weight of 20.5 kDa. Purified by chromatography.
Batch number	E2286-1-1
Batch molecular weight	20.5 kDa

Method	QC requirement	QC Result
HPLC	Reverse Phase HPLC shows >95% purity	Meets specification: >95.0%
Physical appearance	White lyophilised powder solid	Meets specification: White solid
Solubility	Soluble in aqueous buffer	Meets specification

Produced by

Signature

Passed by

Signature

Date

22/06/2023

This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use.